


Product Line up



Product Evaluation Criteria and Acceptable Quality Range


Model Name	MCL-1E	MDL-2E	MVL-3E	UCL-1E	UDL-2E	UDL-3E	UTL-4E	MPC-1E	MPD-2E	MPV-3E
BRAND	Diamond No.1 Meline Classic Meline Plus Classic Admedic Classic No.1 Admedic Plus No.1	Diamond No.2 Meline Deep Meline Plus Deep Admedic Deep No.2 Admedic Plus No.2	Diamond No.3 Meline Volume Meline Plus Volume Admedic Volume No.3 Admedic Plus No.3	Meline Ultra Classic Diamond Plus No.1	Meline Ultra Deep Diamond Plus No.2	Meline Ultra Volume Diamond Plus No.3	Diamond No.4 Meline Ultra Contour Diamond Plus No.4	Meline Pure Classic	Meline Pure Deep	Meline Pure Volume
pH	6.8 ~ 7.2	6.8 ~ 7.2	6.8 ~ 7.2	6.5~7.5	6.5~7.5	6.5~7.5	6.5~7.5	6.5~7.5	6.5~7.5	6.5~7.5
HA content	20mg/mL	20mg/mL	20mg/mL	24mg/mL	24mg/mL	24mg/mL	24mg/mL	20mg/mL	20mg/mL	20mg/mL
Lidocaine content	3.0mg/mL	3.0mg/mL	3.0mg/mL	3.0mg/mL	3.0mg/mL	3.0mg/mL	3.0mg/mL	N/A	N/A	N/A
Residual BDDE	≤2ppm	≤2ppm	≤2ppm	≤2ppm	≤2ppm	≤2ppm	≤2ppm	≤2ppm	≤2ppm	≤2ppm
Osmorality	270~300mmol/kg	270~300mmol/kg	270~300mmol/kg	250~350 mmol/kg	250~350 mmol/kg	250~350 mmol/kg	250~350 mmol/kg	250~350 mmol/kg	250~350 mmol/kg	250~350 mmol/kg
Endotoxin	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL	< 0.5 EU/mL
Viscoelasticity	30 Pa ≤ G' < 80 Pa	80 Pa ≤ G' < 150 Pa	150 Pa ≤ G' < 300 Pa	40 Pa ≤ G' < 90 Pa	100 Pa ≤ G' < 180 Pa	190 Pa ≤ G' < 290 Pa	300 Pa ≤ G' < 360 Pa	40 Pa ≤ G' < 90 Pa	100 Pa ≤ G' < 180 Pa	190 Pa ≤ G' < 290 Pa
Injectability (25G)	N/A	N/A	18 ~ 22 N	N/A	N/A	9 ~ 15 N	10 ~ 16 N	N/A	N/A	8 ~ 14 N
Injectability (27G)	10 ~ 14 N	14 ~ 18 N	10 ~ 14 N	11 ~ 17 N	14 ~ 20 N	19 ~ 25 N	20 ~ 26 N	10 ~ 16 N	13 ~ 19 N	17 ~ 23 N
Injectability (30G)	16 ~ 24 N	N/A	N/A	20 ~ 26 N	N/A	N/A	N/A	19 ~ 25 N	N/A	N/A

Understanding of Wrinkle improvement mechanism




Structure of Younger Skin and Aged Skin

Understanding of Wrinkle improvement mechanism



Before

In general, injected dermal fillers exist in the dermis for 6 to 18 months without chemical changes, and support the lost area to restore volume. However, while they exist in the body, they are gradually broken down into water by the body's metabolic processes and are excreted from the body in the form of moisture.



After

SAFETY - Raw materials

Using **Korea FDA approved materials** such as Non-Animal stabilized Hyaluronic Acid, Lidocaine, and BDDE to ensure quality and safety → **Minimize human allergy incidence**

BLOOMAGE BIOTECHNOLOGY CORP., LTD.
NO.678 TIANCHEN ST., HIGH-TECH DEVELOPMENT ZONE,
JINAN, CHINA 250101

CERTIFICATE OF ANALYSIS 21-FX032204

Product name: **Sodium hyaluronate** Standard: **Ph.Eur.10.0**
Batch No.: 21022801 Manufacture date: Feb.26, 2021
Batch size: 31.28 kg Retest date: Feb.26, 2024
Intrinsic viscosity: **2.2~3.0 ml/kg** Origin: **Bacterial fermentation**

Grade: HA-EP-N 2.5 Intended use: For the manufacture of paracental preparations, including intra-ocular or intra-articular preparations

Items	Acceptance criteria	Results
1. Characters	White or almost white powder or fibrous aggregate	White powder
2. Identification		
A. Infrared absorption	Consistent with the reference substance spectrum of sodium hyaluronate	Complies
B. Reaction of sodium		Positive
3. Appearance of solution	Clear	Clear
A _w ≤ 0.01		0.00
4. pH	5.5~7.0	6.2
5. Intrinsic viscosity	2.2~3.0 ml/kg	2.62 ml/kg
6. Nucleic acids	A _w ≤ 0.1	0.01
7. Protein	≤ 0.1%	< 0.0004%
8. Chlorides	≤ 0.1%	< 0.1%
9. Iron	≤ 20 ppm	< 2 ppm
10. Loss on drying	≤ 10.0%	6.8%
11. Microbial contamination		
TAMC	≤ 100 cfu/g	< 20 cfu/g
TYMC	≤ 100 cfu/g	< 20 cfu/g
12. Bacterial Endotoxins	< 0.04 IU/mg	< 0.04 IU/mg
13. Residual solvents (Ethanol)	≤ 5000 ppm	< 24 ppm
14. Assay	95.0%~105.0% (dried substance)	100.0%

Conclusion: The product complies with the standard of **Ph.Eur.10.0**

Reviewed by: [Signature] Reported by: [Signature]

MAHENDRA CHEMICALS
B-1, 217+1192 G 11 D.C. Estate, Newasa, Ahmedabad, 382 330 Gujarat, INDIA
E-mail: mahendrachemicals1978@gmail.com, Telefax: 091-79-22811390

Hyaluronic Acid

Mahendra Chemicals
B-1, 217+1192 G 11 D.C. Estate, Newasa, Ahmedabad, 382 330 Gujarat, INDIA
E-mail: mahendrachemicals1978@gmail.com, Telefax: 091-79-22811390

Ref. SOP No.: QC 28 **CERTIFICATE OF ANALYSIS**

Product Name : Lidocaine HCl EP

Batch No. : MC /LH /0151/22	A.R. No. : 22FP/0285
Mfg. Date : Dec - 2022	Sample Qty : 20 x 5 gm
Exp. Date : Nov - 2027	Sample Drawn : 02/01/2023
Batch Size : 500.0 Kg	Testing Date : 02/01/2023
Mfg. Lic. No. : G/25/296	Completion Date : 04/01/2023

The result of the test / Analysis as per EP-10.0

Sr. No	TESTS	RESULTS OF ANALYSIS	LIMITS
1.	Characters	A white crystalline powder. Very soluble in water, freely soluble in ethanol	white or almost white crystalline powder. Very soluble in water, freely soluble in ethanol
2.	Identification		
	A Melting Point	74°C to 77°C	Between 74 °C to 79 °C
	B IR	Complies	IR spectrum of sample match with IR spectrum of standard.
	C Reaction of alcoholic Pot. Hydroxide solution	Green color is produced	A green color must be produced
3.	D Reaction of chloride.	Complies	Must comply.
	Appearance of solution	The solution is clear and colourless	Clear and colorless solution
4.	pH	5.11	Between 4.0 to 5.5
5.	Related substance		
	Impurities A	A = Not detected	NMT 0.01 %
6.	Unspecified impurities	Un sp = Not detected	NMT 0.10 %
	Total impurities	Total = Nil	NMT 0.5 %
7.	Sulphated Ash	0.03%	N.M.T. 0.1%
8.	Water (by KF)	6.62%	Between 5.5% and 7.0% w/w
9.	Assay (On anhydrous basis)	99.83%	99.0% To 101.0%
10.	Residual solvent		
	Acetone	Acetone : 2033.5 ppm	NMT 5000 ppm
	Toluene	Toluene : Not detected	NMT 890 ppm

Remarks: This product complies with EP 10.0 Specification.

Analysis By	Checked By	Approved By
Sign: [Signature]	Sign: [Signature]	Sign: [Signature]
Date: 02/01/23	Date: 02/01/23	Date: 02/01/23
Name: A.J. Patel	Name: V.R. Panchal	Name: H.R. Trivedi
Designation: QC Officer	Designation: QC Manager	Designation: QA Manager

Lidocaine

Sigma-Aldrich

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sigma.com
Outside USA: eurtechserv@sigma.com

Certificate of Analysis

Product Name: **1,3-Bis(2-hydroxyethyl)urea - ≥95%**

Product Number: 230892
Batch Number: MKCNT836
BRN#: ALLK4CH
CAS Number: 2425-79-8
MOL Number: MFCC00055146
Formula: C10H18O4
Formula Weight: 202.25 g/mol
Quality Release Date: 17 FEB 2021

Test	Specification	Result
Appearance (Color)	Colorless	Colorless
Appearance (Form)	Liquid	Liquid
Infrared Spectrum	Conforms to Structure	Conforms
Purity (GC)	≥ 95.0 %	97.9 %

Michael Grady, Manager
Quality Control
Milwaukee, WI US

BDDE

Why DIAMOND?



SAFETY - Purification process

Through Separation Refining Technology, remove of residue of cross-linkers by **6 times in 12 days** Dialysis process before sterilization



Enhance safety, reduce side effects

Item	Effect	Diamond	Company A	Company B
Endotoxin(EU/mL)	Inflammation	<0.03	0.5	0.5
BDDE(ppm)	Allergic Reaction	0	<2	<2
Protein Load(ppm)	Allergic Reaction	<0.03	<0.06	<0.05
pH	Pain causing	7.0	7.2	7.3

KTR KOREA TESTING & RESEARCH INSTITUTE

TEST REPORT

38, Cyeplwot-ro, Gaseong-si, Gyeonggi-do, 12670, Korea TEL: 82-(0)-2194-0911 FAX: 82-(0)-2504-1006

Report No : MSK-2023-081343 Receipt Date : 2023.07.05
 Representative : Eyoung Whan Lee Test Completion Date : 2023.07.24
 Company name : Bio Standard Co., Ltd.
 Address : #501, Venture Industry-Academic Cooperation Bldg., 20, Haseo-ro 18beon-gil, Bundang-seo, Seongnam-si, Gyeonggi-do, Republic of Korea
 Sample name : Gynipsothamn, Inactivated or MOC2023

Test Results

TEST ITEM	UNIT	SAMPLE	RESULT	TEST METHOD
TOXICITY - BDE* residual test I (Sublethal, Single dose 0.05g)	mg/kg		0.0 (By the client/Amount)	
-Method (Detection Limit (MDL))				
-BDE: 1 mg/kg				

- Test period : Jul.19.2023 - Jul.24.2023

Approved : By the client - The test methods
 Approved : Form of appearance or shape
 Approved : Lot of lot number

- Usage of Report : For Reference Use

NOTE : 1. The test results of this test report are only limited to the samples and sample names provided by the client and do not guarantee the quality of all products of the client. You can check website (www.ktr.or.kr) or QR code to verify the authenticity of the certificate.
 2. This test report shall be used only within the purpose of its defined usage and shall not be used for public relation, advertisement and so on.
 3. This test report is only valid when printed on KTR original report paper with hologram and when re-issued by KTR. The copy and the electronic file of the test report are only for reference.

Prepared by Sangwon Kim
 Tel : 02-2602-3888

Reviewed by Seung Ho Choi
 Technical Manager
 Tel : 02-2602-3881

2023.07.24

Korea Testing & Research Institute

President *Kim Hyun Cheol*

1 of Total 1 Page(s)

KTR KOREA TESTING & RESEARCH INSTITUTE

QR Code for Integrity

Measurement data from Korea Testing & Research Institute

Why DIAMOND?



SAFETY - Purification process

The BDDE residue test results for BioStandard products conducted by KTR (the Korea Testing & Research Institute) shows 0 ppm, confirming the safety of BioStandard products.

As a Korean government agency, the results of tests conducted by KTR are recognized for their fairness and objectivity.

▶ ▶ ▶
Measurement data
from Korea Testing &
Research Institute

BEYOND ASIAN HUB, TOWARD GLOBAL WORLD

KTR
KOREA TESTING & RESEARCH INSTITUTE

TEST REPORT

98, Goyukwon-ro, Gwacheon-si, Gyeonggi-do, 13810, Korea TEL 82-2-2164-0011 FAX 82-2-2634-1008

Report No : MSK-2023-001343 Receipt Date : 2023.07.05.
Representative : Byung Whan Lee Test Completion Date : 2023.07.24.
Company name : Bio Standard Co.,Ltd.
Address : #501, Venture Industry-Academic Cooperation Bldg., 20, Hoseo-ro 79beon-gil, Baebang-eup, Asan-si, Chungcheongnam-do, Republic of Korea
Sample name : Graft/prosthesis, biomaterial(Lot: MVL2305)

Test Results				
TEST ITEM	UNIT	SAMPLE	RESULT	TEST METHOD
Injection - BDDE residue test (Butanediol diglycidyl ether)	mg/kg		0.0	By the client(Annex#1)

-Method Detection Limit (MDL)-
BDDE: 1 mg/kg
- Test period : Jul.19.2023 ~ Jul.24.2023

Company name : Bio Standard Co.,Ltd.
Address : #501, Venture Industry-Academic Cooperation Bldg., 20, Hoseo-ro 79beon-gil, Baebang-eup, Asan-si, Chungcheongnam-do, Republic of Korea
Sample name : Graft/prosthesis, biomaterial(Lot: MVL2305)

Test Results				
TEST ITEM	UNIT	SAMPLE	RESULT	TEST METHOD
Injection - BDDE residue test (Butanediol diglycidyl ether)	mg/kg		0.0	By the client(Annex#1)

Sanghwoon Kim
Prepared by Sanghwoon Kim
Tel : 02-2092-3898

Seung-Ho Choi
Reviewed by Seung-Ho Choi
Technical Manager
Tel : 1577-0091(ARS ①-④)

2023.07.24.

Korea Testing & Research Institute

President Kim Hyun cheol

1 of Total 1 Page(s)

KTR KOREA TESTING & RESEARCH INSTITUTE

QR Code for forgery


AA1710 Y 2023

Why DIAMOND?




SAFETY - pH & Osmolality

► Comparative analysis of pH



► Comparative analysis of Osmolality




The quality standards closest to the human body's average hydrogen ion concentration **pH 7.0** and osmotic pressure **300 mmol/kg**



- ✓ Minimize human body reactions
- ✓ Reduce side effects

EFFECTIVENESS - Longer Longevity with Shine Technology

SHINE Technology is developed by Bio Standard. This process consists of **3 steps cross-linking** that can last in the body for 6 to 18 months.




Stabilized Hyaluronic acid & Improved Novel Efficacy - SHINE Technology

EFFECTIVENESS - Optimal Viscoelasticity

The optimal viscoelasticity is designed by combining the appropriate viscosity and elasticity for each product type. **Optimal viscoelasticity** by combining the appropriate viscosity and elasticity for each product type







- ✓ **Best for skin generation effect**
- ✓ **Excellent volume effect for a long time**



Comparative analysis of Viscoelasticity between DIAMOND Dermal Fillers

EFFECTIVENES - Optimal viscoelasticity

Injection location	Fine wrinkles around the eyes, earlobes, lips, between the eyes	Lips, nasolabial folds, nose,	Deep nose, nasolabial folds, smile line, forehead	Cheekbones, cheeks, chin, forehead
Model	DIAMOND No.1 MCL-1E	DIAMOND No.2 MDL-2E	DIAMOND No.3 MDL-3E	DIAMOND No.4 UTL-4E
Shape				
Elasticity	Low	Medium	High	Very High
Viscosity	High	Medium	Low	Low
Longevity	6 Months	12 months	15 months	18 months or more

PRODUCT INTRODUCTION – Why DIAMOND?



EFFECTIVENES - Applying Area



Diamond No.1 (MCL-1E)
Fine wrinkles around eyes
Between Eyebrows
Earlobes, Lip



Diamond No.2 (MDL-2E)
Lips, Nasolabial folds,
Nose tip



Diamond No.3 (MVL-3E)
Deep nose, Nasolabial
folds, Smile line,
Forehead



Diamond No.4 (MVL-3E)
Cheek, Chin, Forehead


SATISFACTION - Needle

CE-marked & advanced
ultra-thin wall Needle



Ultra Thin wall


Competitor Thin wall



The difference in the inner diameter of the injection needle due to the difference in the thickness of the injection needle wall

Due to the thinner wall of needle, the outer diameter of needle is the same as that of the competitor's, but the inside diameter is larger, so the **injection pressure is low**. This makes the procedure easier.

SATISFACTION - Even Extrusion force



Compared to competing products, Diamond dermal fillers offer less variation in injection pressure and greater consistency.

Constant injection pressure ensures the practitioner can perform a stable, comfortable, and precise procedure.

Lips

Before



After



Nasolabial folds


Before



After



Nose



Dark circle



INTRODUCTION

Eulji University Medical Center research team in Korea conducted a clinical trial comparing the efficacy of Meline Volume No.3 **DIAMOND No.3**(MHA), with Restylane LYFT (BHA), and it was officially recognized by the **Korean Society of Aesthetic Plastic Surgery**.

ORIGINAL ARTICLE

<https://doi.org/10.14730/aaps.2023.00836>
Arch Aesthetic Plast Surg 2023;29(3):141-146
pISSN: 2184-0831 eISSN: 2288-9537

aaps
Archives of
Aesthetic Plastic Surgery

Comparison of the efficacy and safety between a new monophasic hyaluronic acid filler and a biphasic hyaluronic acid filler in correcting facial wrinkles

Jung Soo Lim¹, Geon Hwi Kim²,
Jong Hun Lee²

¹Vanny Plastic Surgery, Seoul, Korea
²Department of Plastic and Reconstructive Surgery, Nivnon Eulji Medical Center, Eulji University School of Medicine, Seoul, Korea

Background The longstanding and common use of hyaluronic acid (HA) has driven the expanded development of various commercial HA fillers. However, differences in the components of these HA fillers lead to variations in their effect. We compared the in vivo safety and efficacy of biphasic HA (BHA) and a new monophasic HA (MHA) for improving facial wrinkles. We investigated differences in outcomes after their injection into nasolabial folds (NLFs) using the Wrinkle Severity Rating Scale (WSRS), patient satisfaction using the Global Aesthetic Improvement Scale (GAIS), and pain using a visual analog scale (VAS). We also performed a safety assessment of the two fillers.

Methods This matched-pair, double-blind, randomized study compared the degree of temporal wrinkle improvement in the NLFs of 91 participants using the BHA filler versus the new MHA filler. Safety and efficacy were compared at 8 and 24 weeks.

Results At 24 weeks after application, the average WSRS scores were 2.17±0.72 (BHA) and 2.07±0.71 (MHA) (P=0.034). The average GAIS scores, as measured by a treating investigator at 8 weeks and 24 weeks, were 0.94±0.76 (BHA) and 0.98±0.78 (MHA) at 8 weeks (P=0.181), and 0.44±0.64 (BHA) and 0.49±0.69 (MHA) at 24 weeks (P=0.103). The VAS pain score was 0 points at 30 minutes after filler application in both groups.

Conclusions Both the BHA filler and the new MHA filler were safe and effective for improving facial wrinkles in NLFs, but the new MHA filler was more effective for the cosmetic improvement of wrinkle severity than the BHA filler.






Keywords Skin aging / Hyaluronic acid / Dermal fillers / Efficacy / Safety

This work was supported by a grant from the Korea Medical Device Development Fund, funded by the Korean government (the Ministry of Science and ICT, the Ministry of Trade, Industry and Energy, the Ministry of Health and Welfare, the Ministry of Food and Drug Safety) (Project Number: 1711138448, KMDF-PT-20200901_0215).

Source :
<https://doi.org/10.14730/aaps.2023.00836>

Monophasic monodensified fillers such as Meline No. 3 **DIAMOND No.3** with lidocaine (Bio Standard Co., Ltd) share the **characteristics of both biphasic and monophasic polydensified fillers**; therefore, they can exhibit clumping. Overall, however, they spread homogeneously throughout the implanted tissue.

MONOPHASIC VS BIPHASIC

	Monophasic	Biphasic
Structure	 <p> Hyaluronic acid cross link gel</p>	 <p> Hyaluronic acid solution  Hyaluronic acid cross link gel</p>
Product	Juvederm Diamond- Meline Neuramis	Restylane Yvoire Perfectha

MONOPHASIC VS BIPHASIC

	Monophasic	Biphasic
Structure	Homogeneous single-phase product of High cross-linked gel	Two-phase product that mixes cross-linked gel and Non cross-linked gel.
Viscosity	Higher	Lower
Migration	Less Uniform particles penetrate and fill the spaces between collagen and elastin fibers, resulting in very strong cohesion.	More Its larger agglomerate cannot easily penetrate between collagen and elastin fibers, so it is less cohesive.
HA distribution (Spreading)	Homogeneous (Even)	Heterogeneous (Uneven)
Molding	Easier to mold as desired with homogeneous HA distribution.	Hard to mold as desired with uneven heterogeneous HA distribution
Injectability	Easier with even Injection Pressure	Harder with uneven Injection Pressure
Elasticity	Slightly Lower	Higher
Volumization	Slightly less	Superior
Product Classification by Cross-linking	Juvederm DIAMOND_Meline Neuramis	Restylane Yvoire Perfectha

CONCLUSION


Hybrid filler with both monophasic and biphasic properties



Monophasic monodensified fillers such as Meline No. 3_ **DIAMOND No.3** with lidocaine (Bio Standard Co., Ltd) share the **characteristics of both biphasic and monophasic polydensified fillers**; therefore, they can exhibit clumping. Overall, however, they spread homogeneously throughout the implanted tissue.

CONCLUSION

1. DIAMOND No.3 has **lower persistence of pain** compared to Restylane LYFT.




CONCLUSION

Outstanding wrinkle improvement: DIAMOND No.3 showed approximately **0.10 better efficacy** than Restylane LYFT in WSRS (Wrinkle Severity Rating Scale) values.


Table 1. comparison of the efficacy of BHA and MHA fillers for wrinkle improvement

Variable	Evaluator	Time	Mean±SD		Range		P-value
			BHA	MHA	BHA	MHA	
Average WSRS score (n=91)	Independent evaluators	8 wk	2.11±0.72	2.01±0.71	1.0 ~ 0.4		0.055
		24 wk	2.17±0.71	2.07±0.71	1.0 ~ 0.4		0.034



※ WSRS : Wrinkle Severity Rating Scale

Lower WSRS scores indicate superior level of wrinkle improvement



Wrinkle Severity Rating Scale Evaluation

WSRS is evaluated on a total of 5 levels.

The more serious the wrinkle, the higher the score.

The shallower the wrinkle, the lower the score.

In other words, a lower WSRS evaluation score means a better wrinkle condition means that the wrinkles have improved more after injection.

Types of dermal fillers by main ingredient

Main Ingredient	Product	Longevity	Features
Collagen	Zyderm	6 Months	Animal-derived protein extract natural polymer <ul style="list-style-type: none"> - Maintains skin tissue structure - High allergy-causing rate
Hyaluronic acid	Restylene Diamond	6~18 Months	Natural polymer material purified by fermenting microorganisms <ul style="list-style-type: none"> - Wrinkle removal, volume recovery, and skin tone improvement with moisture storage function - Easily can be removed by enzyme (Hyaluronidase), if not satisfied with the results or a medical accident occurs. - More than 80% share due to safety and ease of use
Calcium Hydroxyapatite	Radiesse	18 Months	<ul style="list-style-type: none"> - Cannot be decomposed in the body, so surgical removal is required - A large amount remains, so nodules or limp may occur when a large amount is injected
Poly-L-lactic acid (PLLA)	Scultra	12~24 Months	Promotes collagen production for fine wrinkles <ul style="list-style-type: none"> - At least 4 times injections at 4-week intervals are required - Continuous message is required - High incidence of granulomas with lumpy is found

Types of dermal fillers by main ingredient

Main Ingredient	Product	Longevity	Features
Polymethylmethacrylate (PMMA)	Artefill	24 Months	<ul style="list-style-type: none"> - Induces collagen production as an artificial bone component - Synthetic polymer material used in bovine collagen 80% + PMMA 20%, etc. - Crystals can be created under the skin - Cannot be decomposed in the body, surgical operation required for correction - Absolutely prohibited for use on eyes or forehead, etc.
Poly-caprolactone (PCL)	Elansse	24 Months	<ul style="list-style-type: none"> - Plant-derived synthetic polymer material - Initial volume effect is small, so multiple treatments are required - After the treatment, collagen production is induced - PCL component itself is safe because it is degraded and absorbed, but the collagen produced by PCL can only be corrected and removed through surgery. - Absolutely prohibited for use on the eyes or forehead
Polyacrylamide	Aquamid	24 Months	<ul style="list-style-type: none"> - Particle surface is uneven and granulomas occur, - Severe movement in the body without settlement at the injected parts. - No decomposition in the body, surgical operation required for correction of procedure - Absolutely prohibited for use in the eyes or between the eyebrows

Competitors



Position

	Big Player			Middle Tier	Small Player	Illegal
	Overseas	Korean	Legal Chinese	Biostandard		
Brand	Restylane Juvederm Belotero Perfectha	Neuramis Elavie Yvoir Eptq Revolax(Chaeum)	IMEIK Technology Development / Bloomage Biotechnology / Shanghai Haohai Biological Technology	DIAMOND Sardenya Elasty Replengen Regenovue HyalDew Richesse Rejeunesse	OEM BRAND from Korea, China	A lot from Korea, China, the Third World
Distributor Price	Very High 130 USD ~	High 16 ~ 25USD	Similar to or higher than Korean Products	Low, ~13 USD	~ 7 USD	Very low, less than 5 USD